



Recombinant Human Protein FAM83A (FAM83A)

Product Code	CSB-EP801236HU-B
Abbreviation	FAM83A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86UY5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSRSRHLGKI RKRLEDVKSQ WVRPARADFS DNESARLATD ALLDGGSEAY WRVLSQEGEV DFLSSVEAQY IQAQAREPPC PPDTLGGAEA GPKGLDSSSL QSGTYFPVAS EGSEPALHLS WASAEKPYLK EKSSATVYFQ TVKHNNIRDL VRRICITRTSQ VLVILMDVFT DVEIFCDILE AANKRGVFVC VLLDQGGVKL FQEMCDKVQI SDSHLKNISI RSVEGEIYCA KSGRKFAQI REKFIISDWR FVLSGSYSFT WLCGHVHRNI LSKFTGQAVE LFDEEFRHLY ASSKPVMGLK SPRLVAPVPP GAAPANGRLS SSSGSASDRT SSNPFSGRSA GSHPGTRSVS ASSGPCSPAA PHPPPPPRFQ PHQGPWGAPS PQAHLSPRPH DGPPAAVYSN LGAYRPTRLQ LEQLGLVPRL TPTWRPFLQA SPHF
Source	E.coli
Target Names	FAM83A
Protein Names	Recommended name: Protein FAM83A Alternative name(s): Tumor antigen BJ-TSA-9 Tumor-specific gene expressed in prostate protein
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.