



Recombinant Human Protein FAM124A

(FAM124A)

Product Code	CSB-MP801237HU
Abbreviation	FAM124A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86V42
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDPKAGGGGE EDDCVDSGAE TGGSDYSHLS STSSELSVEE AQDPFLVSIH IADPGESQP LQEAIDNVLA WIHPDLPLFR VSERRASRRR RKPPKGAQPA LAVVLFLQEE YGEEQILQLH RTLQQPPWRH HHTEQVHGRF LPYLPCSQDF FTLAPGTPLW AIRPVHYGKE IVRFTVYCRY DNYADSLRFY QLILRRSPSQ KKADFCIFPI FSNLDVDIQF SLKRLPCDQC PVPTDSSVLE FRVRDIGELV PLLPNPCSPI SEGRWQTEDH DGNKILLQAQ RVHKKFPKPG RVHHASEKKR HSTPLPSTAV PSHTPGSSQQ SPLNSPHPGP IRTGLPPGHQ QEFAGRANST PNPPWSFQRS KSLFCLPTGG PSLASSAEPQ WFSNTGAPGH RASEWRHGHL LSIDDLEGAQ ETDVDTGLRL SSSDLVVSA YSAPSRFCST VETPLPSERC SSHWAAHKDS REGPLPTVSR VTTEASWASL PFFT KRSSSS SATARAAPPA PSTSTLTDSS PQLPCDTPKV KQTDGDMPPP PGSAGPGDND MEEFYI
Source	Mammalian cell
Target Names	FAM124A
Protein Names	Recommended name: Protein FAM124A
Expression Region	1-546
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.