



# Recombinant Human Vacuolar protein-sorting-associated protein 36 (VPS36)

<b>Product Code</b>	CSB-MP801240HU
<b>Abbreviation</b>	VPS36
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q86VN1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDRFVWTSGL LEINETLVIQ QRGVRIYDGE EKIKFDAGTL LLSTHRLIWR DQKNHECCMA ILLSQIVFIE EQAAGIGKSA KIVVHLHPAP PNKEPGPFQS SKNSYIKLSF KEHGQIEFYR RLSEEMTQRR WENMPVSQSL QTNRGPQPGR IRAVGIVGIE RKLEEKRKET DKNISEAFED LSKLMIKAKE MVELSKSIAN KIKDKQGDIT EDETIRFKSY LLSMGIANPV TRETYGSGTQ YHMQLAKQLA GILQVPLEER GGIMSLTEVY CLVNRARGME LLSPEDLVNA CKMLEALKLP LRLRVFDSGV MVIELQSHKE EEMVASALET VSEKGSALTSE EFAKLVGMSV LLAKERLLLA EKMGHLCRDD SVEGLRFYPN LFMTQS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	VPS36
<b>Protein Names</b>	Recommended name: Vacuolar protein-sorting-associated protein 36 Alternative name(s): ELL-associated protein of 45 kDa ESCRT-II complex subunit VPS36
<b>Expression Region</b>	1-386
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.