



Recombinant Human Vacuolar fusion protein MON1 homolog A (MON1A), partial

Product Code	CSB-MP801242HU
Abbreviation	MON1A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86VX9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	DGYKVVFVRRSPLVLVAVARTRQSAQELAQELLYIY YQILSLLTGAQLSHIFQQKQNYDLRRLLSGSERITDNLQLMARDPSFLMGAAR CLPLAA AVRDTVSASLQQARARSLVFSILLARNQLVALVRRKDQFLHPIDLHLLFNLISS SSFRE GEAWTPVCLPKFNAAGFFHAHISYLEPDTDLCLLLVSTDREFFAVSDCRRRF QERLRKR GAHLALREALRTPYYVAQVGIPDLRHFLYKSKSSGLFTSPEIEAPYTSEEEQE RLLGLY QYLHSRAHNASRPLKTIYYTGPNENLLAWVTGAFELYMCYSPLGTKASAVSAI HKLMRWI RKEEDRLFILTPITY
Source	Mammalian cell
Target Names	MON1A
Protein Names	Recommended name: Vacuolar fusion protein MON1 homolog A
Expression Region	205-555
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.