



Recombinant Human RNA-binding protein MEX3D (MEX3D)

Product Code	CSB-MP801257HU
Abbreviation	MEX3D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86XN8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MPSSLGQPDG GGGGGGGGGG VGAAGEDPGP GPAPPPEGAQ EAAPAPRPPP EPDDAAAALR LALDQLSALG LGGAGDTDEE GAAGDGAAAA GGADGGAAPE PVPPDGPEAG APPTLAPAVA PGSLPLLDPN ASPPPPPPPR PSPPDVFAGF APHPAALGPP TLLADQMSVI GSRKKSVMNT ECVVPSSEH VAEIVGRQGC KIKALRAKTN TYIKTPVRGE EPVFIVTGRK EDVEMAKREI LSAAEHFSII RATRSKAGGL PGAAQGPPNL PGQTTIQVRV PYRVVGLVVG PKGATIKRIQ QRTHTYIVTP GRDKEPVFAV TGMPENV DRA REEIEAHITL RTGAFTDAGP DSDFHANGTD VCLDLLGAAA SLWAKTPNQG RRPPTATAGL RGDTALGAPS APEAFYAGSR GGPSVPDPPG ASPYSGSGNG GFAFGAEGPG APVGTAAPDD CDFGFDFDFL ALDLTVPAAA TIWAPFERAA PLPAFSGCST VNGAPGPPAA GARRSSGAGT PRHSPTLPEP GGLRLELPLS RRGAPDPVGA LSWRPPQGPV SFPGGAAFST ATSLPSSPAA AACAPLDGSA SENSRKPPSA SSAPALAREC VVCAEGEVMA ALVPCGHNLF CMDCAVRICG KSEPECPACR TPATQAIHIF S</p>
Source	Mammalian cell
Target Names	MEX3D
Protein Names	Recommended name: RNA-binding protein MEX3D Alternative name(s): RING finger and KH domain-containing protein 1 RING finger protein 193 TINO
Expression Region	1-651
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.