



Recombinant Human Vacuolar protein sorting-associated protein 37D (VPS37D)

Product Code	CSB-EP801259HU
Abbreviation	VPS37D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86XT2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MYRARAARAG PEPGSPGRFG ILSTGQLRDL LQDEPKLDRI VRLSRKFQGL QLEREACLAS NYALAKENLA LRPRLEMGRA ALAIKYQELR EVAENCADKL QRLEESMHRW SPHCALGWLQ AELEEAEQEA EEQMEQLLLG EQSLEAFLPA FQRGRAL AHL RRTQAEKLQE LLRRRERSAQ PAPTSAADPP KSFPAAAVLP TGAARGPPAV PRSLPPLDSR PVPPLKGGSPG CPLGPAPLLS PRPSQPEPPH R
Source	E.coli
Target Names	VPS37D
Protein Names	Recommended name: Vacuolar protein sorting-associated protein 37D Alternative name(s): ESCRT-I complex subunit VPS37D Williams-Beuren syndrome chromosomal region 24 protein
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.