



Recombinant Tomato pseudo-curly top virus Replication enhancer protein (C3, L3)

Product Code	CSB-MP801533TDAC
Abbreviation	C3, L3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q88890
Product Type	Recombinant Protein
Immunogen Species	Tomato pseudo-curly top virus (TPCTV)
Purity	>85% (SDS-PAGE)
Sequence	MDSRTGELLT YTQCQTGCYF DDIKNPIYFR LQNTTTMATG NKIITLQIRA NHNLRKVLGL QICWLNLVVL SVSPRLSGEA LLDRFKFHVM KHLNNLGVIS INNVIRSINH FIDLFDQRVF YAPSFYYNIK MRLY
Source	Mammalian cell
Target Names	C3
Protein Names	Recommended name: Replication enhancer protein Short name= REN Alternative name(s): 16.1 kDa protein Protein C3 Protein L3
Expression Region	1-134
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.