



Recombinant *Oryza sativa* subsp. japonica B3 domain-containing protein Os04g0386900 (Os04g0386900, LOC_Os04g31730)

Product Code	CSB-MP801727OFG
Abbreviation	Os04g0386900, LOC_Os04g31730
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XLP4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MRAATALPSI PSSSSPSPMA SDPTELRCSS PESSGDAGAE DPAAVDAAEE SGGEGGSGHI AAGTEAAPPR PPEPEPEKVA RHGVLPLL GK PYFTCIMCKS HVQPPFQVVV PRSFAPLLPS RTTPATLSWR GRSWGM RFTG GRLIQRLEAG WRGFAVDNDL RLG DGCVFEL LVGGGGGEQER VEFRVQVLRA EIPARIRGRA GGYTSATPIV ID
Source	Mammalian cell
Target Names	Os04g0386900
Protein Names	Recommended name: B3 domain-containing protein Os04g0386900
Expression Region	1-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.