



Recombinant *Oryza sativa* subsp. *japonica* MADS-box transcription factor 17 (MADS17)

Product Code	CSB-MP801740OFG
Abbreviation	MADS17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XUN2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGRGRVELKR IENKINRQVT FSKRRNGLLK KAYELSVLCD AEVALIIFSS RGKLYEFGSA GINKTLEKYN SCCYNAQGSN SALAGGEHQS WYQEMSRLKT KLECLQRSQR HMLGEDLGPL SIKELQQLEK QLEYSLSQAR QRKTQIMMEQ VDDLRRKERQ LGELNKQLKN KLEAEADSSN CRSAIQDSWV HGTVVSGGRV LNAQPPPDID CEPTLQIGYY QFVRPEAANP RSNNGGGGDQN NNFVMGWPL
Source	Mammalian cell
Target Names	MADS17
Protein Names	Recommended name: MADS-box transcription factor 17 Alternative name(s): NMADS3 OsMADS17 RMADS213
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.