



Recombinant *Oryza sativa* subsp. japonica Acidic leucine-rich nuclear phosphoprotein 32-related protein 2 (Os03g0668900, LOC_Os03g46600)

| | |
|--------------------------|---|
| Product Code | CSB-BP801747OFG |
| Abbreviation | Os03g0668900, LOC_Os03g46600 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q7Y180 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. japonica (Rice) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MADAAAPGED DAAWERAIAA AVKNAPFSAP KTLTLDGAVK STTGRLPSPS LLGRYPSLEE LSVAGARLSS LAGLPRLPAL RRLSLPDNRL SGAASLAAVA ESCGATLRHL DLGNNRFADV AELAPLAPHG VESLDLYQCP VTKAKGYRDK VFALIPSLKF LDGMDAEGND CLDSDEEDE EDEGEEGEG EGDEEEEEEG GEEGEGDEDD EEEGDEEDE EEEGEEEADE EDEAGADEED ESKVANGSKG SSGSAQPNKR KRDEDDANG DN |
| Source | Baculovirus |
| Target Names | Os03g0668900 |
| Protein Names | Recommended name: Acidic leucine-rich nuclear phosphoprotein 32-related protein 2 Alternative name(s): ANP32/acidic nuclear phosphoprotein-like protein 2 |
| Expression Region | 1-272 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.