



Recombinant *Oryza sativa* subsp. japonica Acidic leucine-rich nuclear phosphoprotein 32-related protein 2 (Os03g0668900, LOC_Os03g46600)

Product Code	CSB-MP801747OFG
Abbreviation	Os03g0668900, LOC_Os03g46600
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Y180
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MADAAAPGED DAAWERAIAA AVKNAPFSAP KTLTLDGAVK STTGRLPSPS LLGRYPSLEE LSVAGARLSS LAGLPRLPAL RRLSLPDNRL SGAASLAAVA ESCGATLRHL DLGNNRFADV AELAPLAPHG VESLDLYQCP VTKAKGYRDK VFALIPSLKF LDGMDAEGND CLDSDEEDE EDEGEEGEG EGDEEEEEEG GEEGEGDEDD EEEGDEEDE EEEGEEEADE EDEAGADEED ESKVANGSKG SSGSAQPNKR KRDEDDANG DN
Source	Mammalian cell
Target Names	Os03g0668900
Protein Names	Recommended name: Acidic leucine-rich nuclear phosphoprotein 32-related protein 2 Alternative name(s): ANP32/acidic nuclear phosphoprotein-like protein 2
Expression Region	1-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.