



# Recombinant Maturase K (matK)

<b>Product Code</b>	CSB-BP801765PRV
<b>Abbreviation</b>	matK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7YKM5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Proboscidea louisiana (Louisiana Devil's-claw) (Proboscidea louisianica)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEEIQRYLQL ERSQQHDFLY PLIFQEYIYT FAHDHGFSRS IWSKNRGYDN KSSLLIVKRL ITRMYQQNHF LISLNSDNQN PFWARNKNLY SQIISEGFAF IVEIPFSIRL ISCLEGKKIV KSQLNRSIHS IFPFLEXNFS HLNFLVDILI PHPVHVEILD ASSLHLLRFF LNEYCNWNSL ITPKKASSSF SKINQRLFLF LYNSHVCEYE SIFVFLRNQS SHLRSTSSGV LLERIYFYGK IEHLVNVFVT VKDFQANLWL VKEPCMHYIR YQRKSILASK GTSLFMNKWK CYLVTFWQWH FSLWFHPSRI YINQLSNHSL EFLGYLSSVR MNPSVVRSQI LENAFLINNA IKKFDTLVPI IPLIASLAKA KFCNVLGHP SKPVRADLSD SNIIDRFSCI CRNLSHYHSG SSKKKSlyRI KYILRLSCAR TLARKHKSTV RAFLKRLGSE LLEQFLMSEE DVLFXTFXKA SSTLRGVNNS RIWYVDIISI NDLAXHKSKF
<b>Source</b>	Baculovirus
<b>Target Names</b>	matK
<b>Protein Names</b>	Recommended name: Maturase K Alternative name(s): Intron maturase
<b>Expression Region</b>	1-500
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.