



Recombinant Maturase K (matK)

Product Code	CSB-EP801765PRV-B
Abbreviation	matK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7YKM5
Product Type	Recombinant Protein
Immunogen Species	Proboscidea louisiana (Louisiana Devil's-claw) (Proboscidea louisianica)
Purity	≥85% (SDS-PAGE)
Sequence	MEEIQRYLQL ERSQQHDFLY PLIFQEYIYT FAHDHGFSRS IWSKNRGYDN KSSLLIVKRL ITRMYQQNHF LISLNSDNQN PFWARNKNLY SQIISEGFAF IVEIPFSIRL ISCLEGKKIV KSQNLRSIHS IFPFLEXNFS HLNFLVDILI PHPVHVEILD ASSHLLRFF LNEYCNWNSL ITPKKASSSF SKINQRLFLF LYNSHVCEYE SIFVFLRNQS SHLRSTSSGV LLERIYFYGK IEHLVNVFVT VKDFQANLWL VKEPCMHYIR YQRKSILASK GTSLFMNKWK CYLVTFWQWH FSLWFHPSRI YINQLSNHSL EFLGYLSSVR MNPSVVR SQI LENAFLINNA IKKFDTLVPI IPLIASLAKA KFCNVLGHP SKPVRADLSD SNIIDRF GCI CRNL SHYHSG SSKKKS LYRI KYILRLSCAR TLARKHKSTV RAFLKRLGSE LLEQFLMSEE DVLFXTFXKA SSTLRGVNNS RIWYVDIISI NDLAXHKSKF
Source	E.coli
Target Names	matK
Protein Names	Recommended name: Maturase K Alternative name(s): Intron maturase
Expression Region	1-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.