



Recombinant Human Histone H2B type W-T (H2BFWT)

Product Code	CSB-EP801796HU-B
Abbreviation	H2BFWT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z2G1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLRTEVPRLP RSTTAIVWSC HLMATASAMA GPSSETTSEE QLITQEPKEA NSTTSQKQSK QRKRGRHGPR RCHSNCRGDS FATYFRRVLK QVHQGLSLSR EAVSVMDSLV HDILDRIATE AGRLARSTKR QTITAWETRM AVRLLLPQGM GKLAESEGTK AVLRTSLYAI QQQRK
Source	E.coli
Target Names	H2BFWT
Protein Names	Recommended name: Histone H2B type W-T Alternative name(s): H2B histone family member W testis-specific
Expression Region	1-175
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.