



Recombinant Human Growth inhibition and differentiation-related protein 88 (GIDRP88)

Product Code	CSB-EP801814HU-B
Abbreviation	GIDRP88
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z5L2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MQQESERCRV RARRPDMALY VPKARRGAVL LKTGDEEESC GSPNSVVKEK QKESLSQKE VFKDKPEARL LNINPDRKEH NCREEKKSST KL RMDTCLQK TNRVCSKRG T TESKEVLSQG QQQGAPNAGV ITNAPLQRHF KPKKVEC LEV ETTDVTGHER ILLSQACLEI SEAQVPSKPF QNVEFCDFSR HEPDGEAFED KDLEGRIETD TKVLEILYEF PRVFSSVMKP ENMIVPIKLS SDSEIVQQSM QTSDGILNPS SGGITTT SVP GSPDGVFDQT CVD FEVESVG GIANSTGFIL DQKDTDSIPA TMGHISLSES TNDTVSPVMI RECEKNDSTA DELHVKHEPP DTA VLAHETH RDSGFKNVGD ITNKACMMDT TGMSCSDHVT VDSPYVVAVR IADETSINTR SFSKFGMSA DATPLHVAR S GNDTEDFSNP SACSDIYGES ISSHFTESTG KLIESLS DCA SSLPIKKIAG SNYNTFLDSE LSMLNGTKVL SDSAVGIDLG STGDTTEALH ELRTAEFFKT EEQDDSGSIE FGVSFPDRES SSMETSIEPK ATETSHTEGI TAIEESWESM FNDDGDCLDP RLLQEGILMH IKPENHCSEL SGNTKSRESI QEPRSDYYNH EVPDIDLSDC EFPHVIEIYD FPQEFHTEDL LRVFCSYQKK GFDIKWVDDT HALGVFSSPI TARDALGIKH TMVKIRPLSQ ATRAAKAKAR AYA EFLQPAK ERPETSAALA RRLVISALGV RSKQSKTERE AELKKLQEAR ERKRLEAKQR EDIWEGRDQS TV </p>
Source	E.coli
Target Names	R3HCC1L
Protein Names	Recommended name: Growth inhibition and differentiation-related protein 88 Alternative name(s): Putative mitochondrial space protein 32.1
Expression Region	1-792
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.