



Recombinant Human Coiled-coil domain-containing protein 91 (CCDC91)

Product Code	CSB-EP801816HU-B
Abbreviation	CCDC91
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z6B0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDDDDFGGFE AAETFDGGSG ETQTTSPAIP WAAFPVAVSGV HLSPSSPEIV LDRDHSSSIG CLSSDAIISS PENTHAANSI VSQTIPKAQI QQSTHTHLDI SLFPLGLTDE KSNGTIALVD DSEDPGANVS NIQLQQKISS LEIKLKVSEE EKQRIKQDVE SLMEKHNVLE KGFLKEKEQE AISFQDRYKE LQEKHKQELE DMRKAGHEAL SIIVDEYKAL LQSSVKQQVE AIEKQYISAI EKQAHKCEEL LNAQHQRLLLE MLDTEKELLK EKIKEALIQQ SQEQKEILEK CLEERQRNK EALVSAAKLE KEAVKDAVLK VVEEERKNLE KAHAEERELW KTEHAKDQEK VSQEIQKAIQ EQRKISQETV KAAIIIEQKR SEKAVEEAVK RTRDELIEYI KEQKRLDQVI RQRSLSSLEL FLSCAQKQLS ALIATEPVDI E
Source	E.coli
Target Names	CCDC91
Protein Names	Recommended name: Coiled-coil domain-containing protein 91 Alternative name(s): GGA-binding partner p56 accessory protein
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.