



Recombinant *Xenopus laevis* Eukaryotic translation initiation factor 2 subunit 1 (eif2s1)

Product Code	CSB-YP801857XBE
Abbreviation	eif2s1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZTK4
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	PGLNCRFYQ HKFPEVDDVV MVNVRSAEM GAYVSLLEYN NIEGMILLSE LSRRRIRSIN KLIRIGRNEC VVVIRVDKDK GYIDLSKRRV SPEALKCED KFTKSKTVYS ILRHVAEVLD YTKDEQLDSL FQRTAWVFDE KYKKPGYGAY DAFKNAVSDP DILDGLDLSE DERRVLIDNI NRRLTPQAVK IRADIEVACY GYEGIDAVKD ALRAGLSCST ENMPIKINLI APPRYVMTTT TLERTEGLSV LNQAMSVIKE RIEEKRGVFN VQMEPKVVTD TDETELARQL ERLEKENAEV DGDDDADEME AKTED
Source	Yeast
Target Names	eif2s1
Protein Names	Recommended name: Eukaryotic translation initiation factor 2 subunit 1 Alternative name(s): Eukaryotic translation initiation factor 2 subunit alpha Short name= eIF-2-alpha Short name= eIF-2A Short name= eIF-2alpha
Expression Region	2-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.