



Recombinant Danio rerio Abhydrolase domain-containing protein FAM108C1 (fam108c1)

Product Code	CSB-BP801871DIL
Abbreviation	fam108c1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZVZ7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MPEQGPRMNS FSLGELCWLF CCPPCPSRIA AKLAFLPPEP TYSVHTDPSG ATSLHLTERA DWQYSQRELD AVEVLVTRTS RGNRVGCMFV RCAPASRYTL LFSHGNAVDL GQMCSFYIGL GSRINCNVFS YDYSGYGVST GKPSEKNLYA DIEAAWQVLR NKYGVTPEI ILYGQSIGTV PTVDLASRYE CAAVILHSPL MSGLRVAFPD TRKTYCFDAF PSIDKVSVA SPVLVIHGTE DEVIDFSHGL AIYERCPRAV EPLWVEGAGH NDIELYAQYL ERLKQFITFE LATS
Source	Baculovirus
Target Names	abhd17c
Protein Names	Recommended name: Abhydrolase domain-containing protein FAM108C1 EC=3.-.-
Expression Region	1-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.