



Recombinant Protein sprT (sprT)

Product Code	CSB-MP802803SZB
Abbreviation	sprT
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q83JS6
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MKTSRLPIAI QQAVMRRLLRE KLAQANLKLK RNYPEPKLSY TQRGTSAGTA WLESYEIRLN PVLLLENSA FIEEVVPHL AHLLVWKHFG RVAPHGKEWK WMMESVLGVP ARRTHQFELQ SVRRNTFPYR CKCQEHQLTV RRHNRVVRGE AIYRCVHCGE QLVAK
Source	Mammalian cell
Target Names	sprT
Protein Names	Recommended name: Protein sprT
Expression Region	1-165
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>