



Recombinant 3-phenylpropionate/cinnamic acid dioxygenase subunit alpha (hcaE)

Product Code	CSB-MP802806SZB
Abbreviation	hcaE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q83K39
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MTTPSDLNIY QLIDTQNGRV TPRIYTPDI YQLELERIFG RCWLFLAHES QIPKPGDFFN TYMGEDAVVV VRQKDGSIKA FLNQCRHRAM RVSYADCGNS RAFTCPYHGW SYGINGELID VPLEPRAYPQ GLCKSHWGLN EVPCVESYKG LIFGNWDTSA PGLHDYLGDI AWYLDGMLDR REGGTEIVGG VQKWVINCNW KFPAEQFASD QYHALFSHAS AVQVLGAKDD GSDKRLGDGQ TARPVWETAK DALQFGQDGH GSGFFFTEKP DANVWVDGAV SSYYRETYAE AEQRLGEVRA LRLAGHNNIF PTLSWLNGTA TLRVWHPRGP DQVEVWAFCI TDKAASDEVK AAFENSATRA FGPAGFLEQD DSENVCEIQK LLKGHRARNS KLCLEMGLGQ EKRRDDGIPG ITNYIFSETA ARGMYQRWAD LLSSESWQEV LDKTAAYQQE VMK
Source	Mammalian cell
Target Names	hcaE
Protein Names	Recommended name: 3-phenylpropionate/cinnamic acid dioxygenase subunit alpha EC= 1.14.12.19 Alternative name(s): Digoxigenin subunit alpha
Expression Region	1-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.