



Recombinant Homoserine kinase (thrB)

Product Code	CSB-EP802830SZB
Abbreviation	thrB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q83MH9
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	≥85% (SDS-PAGE)
Sequence	MVKVYAPASS ANMSVGFDVL GAAVAPVDGA LLGDVVTVEA AETFSLNNLG RFADKLPSEP RENIVYQCWE RFCQELGKQI PVAMTLEKNM PIGSGLGSSA CSVVAALMAM NEHCGKPLND TRLLALMGEL EGRISGSIHY DNVAPCFLGG MQLMIEENDI ISQQVPGFDE WLWVLAYPGI KVSTAEARAI LPAQYRRQDC IAHGRHLAGF IHACYSRQPE LAAKLMKDVI AEPYRERLLP GFRQARQAVA EIGAVASGIS GSGPTLFALC DKPDTAQRVA DWLGKNYLQN QEGFVHICRL DTAGARVLEN
Source	E.coli
Target Names	thrB
Protein Names	Recommended name: Homoserine kinase Short name= HK Short name= HSK EC= 2.7.1.39
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.