



Recombinant Arabidopsis thaliana Ribosomal RNA-processing protein 8 (At5g40530)

Product Code	CSB-EP802900DOA-B
Abbreviation	At5g40530
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84JC0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MTTEENKTSR NRKRKRQRNP KPSKEEPIET TPKNQNEKKN QRDTKNQQHG GSSAPSKRPK PSNFLDALRE RLSGGQFRML NEKLYTCSGK EALDYFKEDP QMFDMYHTGY QQQMSNWPEL PVNSIINWLL SNSSSLVVAD FGCGDARIAK SVKNKVFSFD LVSKNPSVIA CDMSNTSLES SSVDVAVFCL SLMGTNYSSY IKEAHRVLRP SGMLLIAEVK SRFDPNNGGA DPKDFVKAVC DLGFTSVLKD FSNKMFILFH FKKKEQMNSN QKIIKWPELK ACLYKRR
Source	E.coli
Target Names	At5g40530
Protein Names	Recommended name: Ribosomal RNA-processing protein 8 EC= 2.1.1.-
Expression Region	1-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.