



# Recombinant Arabidopsis thaliana Uncharacterized protein At4g37920, chloroplastic (At4g37920)

<b>Product Code</b>	CSB-MP802934DOA
<b>Abbreviation</b>	At4g37920
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q84WN0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AEVKSSVE DPMEVEVAEG YTMAQFCDKI IDLFLNEKPK VKQWKTYLVL RDEWNKYSVN FYKRCRIRAD TETDPILKQK LVSLESKVKK IDKEMEKHND LLKEIQENPT DINAI AAKRR RDFTGGEFFRY VTLLSETLDG LEDRDAVARL ATRCLSAVSA YDNTLESVET LDTAQAKFED ILNSPSVDSA CEKIRSLAKA KELDSSLILL INSAYAAAKE SQTVTNEAKD IMYHLYKATK SSLRSITPKE IKLLKYLLNI TDPEERFSAL ATAFSPGDDH EAKDPKALYT TPKELHKWIK IMLDAYHLNK EETDIKEAKQ MSQPIVIQRL FILKDTIEDE YLDKKTIVAD ETPKKEEEDT TIEDFLN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At4g37920
<b>Protein Names</b>	Recommended name: Uncharacterized protein At4g37920, chloroplastic
<b>Expression Region</b>	63-427
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.