



Recombinant C-type lectin domain-containing protein 88 (clec-88)

Product Code	CSB-EP803099CXY-B
Abbreviation	clec-88
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86NG3
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	VT NDIEDASGET PGIVSQITEE QPHQRQKLYN WDYKDLGTIA FEDIPFPTLQ PSQTIDQSEN CPEGWIRYSD SCYWVETELL GFAKAERKCS EKQSTLFVAN SIDEWEAIRG HSKDSFSSWI GLVRFSHYER SEQLPRWQTE GAVNPSKINW LIKPYSPVLN GWSQLANCAA SYKSPASLET ASYTYFYPCT YLFYSICERN STIANSLH
Source	E.coli
Target Names	clec-88
Protein Names	Recommended name: C-type lectin domain-containing protein 88 Alternative name(s): Chondroitin proteoglycan 6
Expression Region	19-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.