



# Recombinant C-type lectin domain-containing protein 88 (clec-88)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP803099CXY   |
| <b>Abbreviation</b>      | clec-88   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q86NG3  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Caenorhabditis elegans  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | VT NDIEDASGET PGIVSQITEE QPHQRQKLYN WDYKDLGTIA FEDIPFPTLQ<br>PSQTIDQSEN CPEGWIRYSD SCYWVETELL GFAKAERKCS EKQSTLFVAN<br>SIDEWEAIRG HSKDSFSSWI GLVRFSHYER SEQLPRWQTE GAVNPSKINW<br>LIKPYSPVLN GWSQLANCAA SYKSPASLET ASYTYFYPCT YLFYSICERN<br>STIANSLH   |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | clec-88   |
| <b>Protein Names</b>     | Recommended name: C-type lectin domain-containing protein 88<br>Alternative name(s): Chondroitin proteoglycan 6   |
| <b>Expression Region</b> | 19-228  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |