



Recombinant Human Ankyrin repeat domain-containing protein 23 (ANKRD23)

Product Code	CSB-BP803108HU
Abbreviation	ANKRD23
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86SG2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDFISIQQLV SGERVEGKVL GFGHGVDPDG AWPSDWRRGP QEAVAREKLK LEEEKKKKLE RFNSTRFNLD NLADLENLVQ RRKKRLRHRV PPRKPEPLVK PQSQAQVEPV GLEMFLKAAA ENQEYLIDKY LTDGGDPNAH DKLHRTALHW ACLKGHSQLV NKLLVAGATV DARDLLDRTP VFWACRGGHL VILKQLLNQG ARVNARDKIG STPLHVAVRT RHPDCLEHLI ECGAHLNAQD KEGDTALHEA VRHGSYKAMK LLLLYGAELG VRNAASVTPV QLARDWQRGI REALQAHVAH PRTRC
Source	Baculovirus
Target Names	ANKRD23
Protein Names	Recommended name: Ankyrin repeat domain-containing protein 23 Alternative name(s): Diabetes-related ankyrin repeat protein Muscle ankyrin repeat protein 3
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is a member of the muscle ankyrin repeat protein (MARP) family and encodes a protein with four tandem ankyrin-like repeats. The protein is localized to the nucleus, functioning as a transcriptional regulator. Expression of this protein is induced during recovery following starvation.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.