



Recombinant Human Zinc finger and BTB domain-containing protein 46 (ZBTB46)

Product Code	CSB-EP803124HU-B
Abbreviation	ZBTB46
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86UZ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNNRKEDMEI TSHYRHLLRE LNEQRQHGV L CDVCVVVEGK VFKAHKNVLL GSSRYFKTLY CQVQKTSEQA TVTHLDIVTA QGFKAIIDFM YSAHLALTSR NVIEVMSAAS FLQMTDIVQA CHDFIKAALD ISIKSDASDE LAEFEIGASS SSSTEALISA VMAGRSISPW LARRTSPANS SGDSAIASCH DGGSSYGKED QEPKADGPDD VSSQPLWPGD VGYGPLRIKE EQVSPSQYGG SELPSAKDGA VQNSFSEQSA GDAWQPTGRR KNRKNKETVR HITQQVEDDS RASSPVPSFL PTSGWPFSSR DSNADLSVTE ASSSDSRGER AELYAQVEEG LLGGEASYLG PPLTPEKDDA LHQATAVANL RAALMSKNSL LSLKADV LGD DGSLLFEYLP RGAHSLSLNE FTVIRKKFKC PYCSFSAMHQ CILKRHMRSH TGERPYPCEI CGKKFTRREH MKRHTLVHSK DKKYVCKVCS RVFMSAASVG IRHGSRRHGV CTDCAGRGMA GPLDHGGGGG EGSPEALFPG DGPYLEDPED PRGEAEELGE DDEGLAPEDA LLADDKDEED SPRPRSPPGG PDKDFAWLS
Source	E.coli
Target Names	ZBTB46
Protein Names	Recommended name: Zinc finger and BTB domain-containing protein 46 Alternative name(s): BTB-ZF protein expressed in effector lymphocytes Short name= BZEL BTB/POZ domain-containing protein 4 Zinc finger protein 340
Expression Region	1-589
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.