



Recombinant Human PNMA-like protein 1 (PNMAL1)

Product Code	CSB-MP803126HU
Abbreviation	PNMAL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86V59
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSKTMAMNLL EDWCRGMEVD IHRSLLVGTI PEDCGQAEIE ETLNGVLSPL GPYRVLNKIF VREENVKAAL IEVGEGVNLS TIPREFPGRG GVWRVVCRRDP TQDAEFLKNL NEFLDAEGRT WEDVVRLQL NHPTLSQNLQH QPPENWAEAL GVLLGAVVQI IFCMDAEIRS REEARAQEAA EFEEMAAWAL AAGRKVKKEP GLAAEVGSAL KAETPNWNNA TEDQHEPTKP LVRRAGAKSR SRRKKQKKNS RQEAVPWKKP KGINSNSTAN LEDPEVGDAE SMAISEPIKG SRKPCVNKEE LALKKPMAC AWKGPREPPQ DARAEAESPG GASESDQDGG HESPPKKKAV AWVSAKNPAP MRKKKKVSLG PVSYYLVDSE DGRKKPVMPK KPGSRREAS DQKAPRGQQP AEATASTSRG PKAKPEGSPR RATNESRKV
Source	Mammalian cell
Target Names	PNMA8A
Protein Names	Recommended name: PNMA-like protein 1
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.