



Recombinant Human Ran-binding protein 3-like (RANBP3L)

Product Code	CSB-EP803130HU-B
Abbreviation	RANBP3L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86VV4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTTIPRKGSS HLPGLSHTCK LKLQEDRRQQ EKSVIAQPIF VFEKGEQTFK RPAEDTLYEA AEPECNGFPT KRVRSSSFTF HITDSQSQGV RKNNVFMTSA LVQSSVDIKS AEQGPVKHSHK HVIRPAILQL PQARSCAKVR KTFGHKALES CKTKEKTNNK ISEGNSYLLS ENLSRARISV QLSTNQDFLG ATSVGCQPNE DKCSFKSCSS NFVFGENMVE RVLGTQKLTQ PQLENDSYAK EKPFKSIPKF PVNFLSSRTD SIKNTSLIES AAAFSSQPSR KCLLEKIDVI TGEETEHNVL KINCKLFIFN KTTQSWIERG RGTLRLNDDTA STDCGTLQSR LIMRNQGSRL LILNSKLWAQ MKIQRANHKN VRITATDLED YSIKIFLIQA SAQDTAYLYA AIHHRLVALQ SFNKQRDVNQ AESLSETAQQ LNCESCDENE DDFIQVTKNG SDPSSWTHRQ SVACS
Source	E.coli
Target Names	RANBP3L
Protein Names	Recommended name: Ran-binding protein 3-like
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.