



Recombinant Human Tetratricopeptide repeat protein 30A (TTC30A)

Product Code	CSB-MP803138HU
Abbreviation	TTC30A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86WT1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAGLSGAQIP DGEFTALVYR LIRDARYAEA VQLLGRELQR SPRSRAGLSL LGYCYYRLQE FALAAECYEQ LGQLHPELEQ YRLYQAQALY KACLYPEATR VAFLLLDNPA YHSRVLRLQA AIKYSEGDLP GSRSLVEQLL SGEGGEESGG DNETDGQVNL GCLLYKEGQY EAACSKFSAT LQASGYQPD LSYNLALAYYS SRQYASALKH IAEIIRGIR QHPELGVGMT TEGFDVRSVG NTLVLHQ TAL VEAFNLKAAI EYQLRNYEVA QETLTDMPPR AEEELDPVTL HNQALMNM DA RPTGFEKLQ FLLQQNPFPP ETEFGNLLLLY CKYEYFDLAA DVLAENAH LT YKFLTPYLYD FLDALITCQT APEEAFIKLD GLAGMLTEQL RRLTKQVQEA RHNRDDEAIK KAVNEYDETM EKYIPVLMAQ AKIYWNLENY PMVEKVFRKS VEFCNDHDVW KLNVAHVLFM QENKYKEAIG FYEPIVKKHY DNILNVSAIV LANLCVSYIM TSQNEEA EEL MRKIEKEEEEQ LSYDDPNRKM YHLCIVNLVI GTLYCAKGN Y EFGISRVIKS LEPYNKLG TDTWYAKRCF LSLLNMSKH MIVIHDSVIQ ECVQFLGHCE LYGTNIPAVI EQPLEEERMH VGKNTVTDES RQLKALIYEI IGW NK
Source	Mammalian cell
Target Names	TTC30A
Protein Names	Recommended name: Tetratricopeptide repeat protein 30A Short name= TPR repeat protein 30A
Expression Region	1-665
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.