



Recombinant Human HORMA domain-containing protein 1 (HORMAD1)

Product Code	CSB-EP803140HU-B
Abbreviation	HORMAD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86X24
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MATAQLQRTP MSALVFPNKI STEHQSLVLV KRLLAVSVSC ITYLRGIFPE CAYGTRYLDD LCVKILREDK NCPGSTQLVK WMLGCYDALQ KKYLRMVVLA VYTNPEDPQT ISECYQFKFK YTNGPLMDF ISKNQSNES MLSTDTKKAS ILLIRKIYIL MQNLGPLPND VCLTMKLFYY DEVTPPDYQP PGFKDGDCEG VIFEGEPMYL NVGEVSTPFH IFKVKVTTTER ERMENIDSTI LSPKQIKTPF QKILRDKDVE DEQEHYTSDD LDIETKMEEQ EKNPASSELE EPSLVCEEDE IMRSKESPD L SISHSQVEQL VNKTSELDMS ESKTRSGKVF QNKMANGNQP VKSSKENRKR SQHESGRIVL HHFDSSSQES VPKRRKFSEPE KEHI
Source	E.coli
Target Names	HORMAD1
Protein Names	Recommended name: HORMA domain-containing protein 1 Alternative name(s): Cancer/testis antigen 46 Short name= CT46 Newborn ovary HORMA protein
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.