



Recombinant Human Sorting nexin-32 (SNX32)

Product Code	CSB-YP803143HU
Abbreviation	SNX32
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q86XE0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	METYAEVGKE GKPCASVDL QGDSSLQVEI SDAVSRDKV KFTVQTKSCL PHFAQTEFSV VRQHEEFIWL HDAYVENEY AGLIIPPAPP RPDFEASREK LQKLGEGDSS VTREEFAKMK QELEAEYLAI FKKTVMHEV FLQRLAAHPT LRRDHNFFVF LEYGQDLSVR GKNRKELGG FLRNIVKSAD EALITGMSGL KEVDDFFEHE RTFLLEYHTR IRDACLRADR VMRAHKCLAD DYIPISAALS SLGTQEVNQL RTSFLKLAEL FERLRKLEGR VASDEDLKLS DMLRYMYRDS QAAKDLLYRR LRALADYENA NKALDKARTR NREVRPAESH QQLCCQRFER LSDSAKQELM DFKSRRVSSF RKNLIELAEL ELKHAKASTL ILRNTLVALK GEP
Source	Yeast
Target Names	SNX32
Protein Names	Recommended name: Sorting nexin-32 Alternative name(s): Sorting nexin-6B
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>