



Recombinant Human Rab11 family-interacting protein 4 (RAB11FIP4)

Product Code	CSB-MP803157HU
Abbreviation	RAB11FIP4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86YS3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAGGAGWSGA PAALLRSVRR LREVFVCGR DPDGFLRVER VAALGLRFGQ GEEVEKLVKY LDPNDLGRIN FKDFCRGVFA MKGCEELLKD VLSVESAGTL PCAPEIPDCV EQGSEVTGPT FADGELIPRE PGFFPEDEEE AMTLAPPEGP QELYTDSPME STQSLEGSVG SPAEKDGGGLG GLFLPEDKSL VHTPSMTTSD LSTHSTTSLI SNEEQFEDYG EGDDVDCAPS SPCPDDETR NVYSDLGSSV SSSAGQTPRK MRHVYNSELL DVYCSQCCKK INLLNDLEAR LKNLKANSPN RKISSTAFGR QLMHSSNFSS SNGSTEDLFR DSIDSCDNDI TEKVSFLEKK VTELENSLT NGDLKSKLKQ ENTQLVHRVH ELEEMVKDQE TTAEQALEEE ARRHREAYGK LEREKATEVE LLNARVQQLE EENTELRTTV TRLKSQTEKL DEERQRMSDR LEDTSLRLKD EMDLYKRMMMD KLRQNRLEFQ KEREATQELI EDLRKELEHL QMYKLDLDCERP GRGRSASSGL GEFNARAREV ELEHEVKRLK QENYKLRDQN DDLNGQILSL SLYEAKNLFA AQTKAQSLAA EIDTASRDEL MEALKEQEEI NFRLRQYMDK IILAILDHNP SILEIKH
Source	Mammalian cell
Target Names	RAB11FIP4
Protein Names	Recommended name: Rab11 family-interacting protein 4 Short name= FIP4-Rab11 Short name= Rab11-FIP4 Alternative name(s): Arfophilin-2
Expression Region	1-637
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.