



# Recombinant Pseudomonas syringae pv. tomato Enolase 1 (eno1)

<b>Product Code</b>	CSB-EP803408FGP-B
<b>Abbreviation</b>	eno1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q886M3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pseudomonas syringae pv. tomato (strain ATCC BAA-871 / DC3000)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAKIVDIKGR EVLDSRGNPT VEADVLLDNG IIGSACAPSG ASTGSREALE LRDGDKSRYM GKGVLKAVAN INGPIRDLLL GKDPVDQKAL DHAMIALDAT ENKASLGANA ILAVSLAAAK AAAQDQDLPL YAHIANLNGT PGVYSMPVPM MNIINGGEHA DNNIDIQEFM IQPVGAKSFA EALRWGTEIF HHLKAVLKAR GLNTAVGDEG GFAPNLASNK EALEAIAEAV ANAGYTLGTD VTLALDCAAS EFYKNGKYTL SEEGEYSSAE FAEYLAELTR KHPIISIEDG LDESDDWDGWK VLTEKIGEKV QLVGDDLFTV NTKILKEGID KNIANSILIK FNQIGTLTET LEAIQMAKAA GYTAISHRS GETEDSTIAD LAVGTSAGQI KTGSLCRSDR VSKYNQLLRI EEQLGSKAVY RGRAEFRG
<b>Source</b>	E.coli
<b>Target Names</b>	eno1
<b>Protein Names</b>	Recommended name: Enolase 1 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 1 2-phosphoglycerate dehydratase 1
<b>Expression Region</b>	1-428
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.