



Recombinant Walleye dermal sarcoma virus Protein ORF-C (orfC)

Product Code	CSB-YP803424WBH
Abbreviation	orfC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q88936
Product Type	Recombinant Protein
Immunogen Species	Walleye dermal sarcoma virus (WDSV)
Purity	>85% (SDS-PAGE)
Sequence	MAWYHQHRWH LDYSIPRQNL QAFLTTITFI DPQFKIQENG LTEGEYKTQI VKQIIPQLCR IPNQNSPPPI WVQGPRIKGD PTWLKINAKF ITLIPKQKG TKNISTKTYL SRLFVIWLQN
Source	Yeast
Target Names	orfC
Protein Names	Recommended name: Protein ORF-C
Expression Region	1-120
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.