



# Recombinant Mouse Complement component C8 beta chain (C8b)

<b>Product Code</b>	CSB-EP803869MO-B
<b>Abbreviation</b>	C8b
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8BH35
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>SVAEAPA PIDCELSTWS SWTACDPCQK KRYRHTYLLR PSQFYGELCD  LSDKEVEDCV TNQPCRSQVR CEGFVCAQTG RCVNRRLLCN  GDNDCGDQSD EANCRRYKN CQREMEQYWA IDRLASGINL FTNTFEGPVL  DHRYYAGGCS PHYILD TNFR KPYNVESYTP QTKCEYEFTL TEYESYSDFE  RLVIEK KTHM FNFTSGFKVD GVMDLGIKVE SNEGKNYVTR TKRFAHTQSK  FLHARSVLEV AHYKLKSRSL MLHYEFLQRV KSLPLEYSYG EYRDLLRDFG  THFITEAVLG GIYEYTLIMN KDAMEQGDYT LSHVTACAGG SFGIGGMVYK  VYVKVGVSAK KCSDIMKEIN ERNKRSTMVE DLVVLVRGGT SEDITALAYK  ELPTPELMEA WGDVAVKYNPA IIKIKAEPLY ELVTATDFAY SSTVKQNLKK  ALEEFQSEVS SCRCAPCRGN GVPVLKGSRC ECICPGGFQG TACEVTYRKD  IPIDGKWSCW SDWSACSGGH KTRHRQCNNP APHKGGSPCS GPASETLCN</p>
<b>Source</b>	E.coli
<b>Target Names</b>	C8b
<b>Protein Names</b>	Recommended name: Complement component C8 beta chain Alternative name(s): Complement component 8 subunit beta
<b>Expression Region</b>	54-589
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	C8 beta is one of the three subunits that comprise the component 8 (C8) of the complement system. C8 participates in the formation of Membrane Attack Complex that results in the lysis of cells. Patients with C8B deficiency are prone to bacteria infection.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.