



# Recombinant Mouse Polycomb group RING finger protein 3 (PcgF3)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP803934MO-B                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Abbreviation</b>      | PcgF3                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |
| <b>Uniprot No.</b>       | Q8BTQ0                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Mus musculus (Mouse)                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | MLTRKIKLWD INAHITCRLC SGYLIDATTV TECLHTFCRS CLVKYLEENN<br>TCPTCRIVIH QSHPLQYIGH DRTMQDIVYK LVPGLQEAEM RKQREFYHKL<br>GMEVPGDIKG EACSAKQHLD PRNGETKADD NSNKETAEEK QEEDNDYHRS<br>DEQVSICLEC NSSKLRLGLKR KWIRCSAQAT VLHLKKFIAK KLNLSFNL<br>DILCNEEILG KDHTLKFVVV TRWRFKKAPL LLHYRPKMDL L                                                                                                                        |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Target Names</b>      | PcgF3                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Names</b>     | Recommended name: Polycomb group RING finger protein 3 Alternative name(s): RING finger protein 3A                                                                                                                                                                                                                                                                                                          |
| <b>Expression Region</b> | 1-241                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Target Details</b>    | This protein contains a C3HC4 type RING finger, which is a motif known to be involved in protein-protein interactions. The specific function of this protein has not yet been determined.                                                                                                                                                                                                                   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |