



# Recombinant Mouse Sialidase-4 (Neu4)

<b>Product Code</b>	CSB-MP803970MO
<b>Abbreviation</b>	Neu4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8BZL1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGPTRVPRRT VLFQRERTGL TYRVPALLCV PPRPTLLAFA EQRLSPDDSH AHRLVLRRTG LTRGSVRWGT LSVLETAVLE EHRSMNPCPV LDEHSGTIFL FFIAVLGHTP EAVQIATGKN AARLCCVTSC DAGLTWGSVR DLTEEAIGAA LQDWATFAVG PGHGVQLRSG RLLVPAYTYH VDRRECFGKI CWTSPHSLAF YSDDHGISWH CGGLVPLNRS GECQLAAVDG DFLYCNARSP LGNRVQALSA DEGTSFLPGE LVPTLAETAR GCQGSIVGFL APPSIEPQDD RWTGSPRNT HSPCFNLRVQ ESSGEGARGL LERWMPRLPL CYPQSRSPEN HGLEPGSDGD KTSWTPECPM SSSDMLQSPT WLLYSHPAGR RARLHMGIYL SRSPLDPHSW TEPWVIYEGP SGYSDLAFLG PMPGASLVFA CLFESGTRTS YEDISFCLFS LADVLENVPT GLEMLSLRDK AQGHCWPS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Neu4
<b>Protein Names</b>	Recommended name: Sialidase-4 EC= 3.2.1.18 Alternative name(s): N-acetyl-alpha-neuraminidase 4 Neuraminidase 4
<b>Expression Region</b>	1-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein belongs to a family of glycohydrolytic enzymes, which remove terminal sialic acid residues from various sialo derivatives, such as glycoproteins, glycolipids, oligosaccharides, and gangliosides. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .