



Recombinant Chicken Endophilin-A3 (SH3GL3)

Product Code	CSB-MP804294CH
Abbreviation	SH3GL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8AXU9
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MSVAGLKKQF HKASQLFSEK ISGAEGTKLD EEFQEMERKI DVTNKAVAE LSKSTEYLQP NPAYRAKLG M LNTMSKIRGQ VKTTGYPQTE GLLGDCMIRY GRELGDSDMF GLALLDAGES MKQMAEVKDS LDINVKQNF I DPLQLLQDKD LKEIGHHLKK LEGRRLDYDY KKKRLGKIPD EEVKQAVEKF EESKELAERS MFNFLNDVE QVSQLAVFVE AALDYHKQST EILEDLQSKL QNRINVASSR PKREFKPKPV ITTTLETGDN QQHNGIAYSS SIKSSGSSMH VDQPCCQALY DFEPENEGEL GFKEGDIITL TNQIDENWYE GMLNGESGFF PHNYVEVMVP LPQ
Source	Mammalian cell
Target Names	SH3GL3
Protein Names	Recommended name: Endophilin-A3 Alternative name(s): Endophilin-3 SH3 domain-containing GRB2-like protein 3 SH3p13
Expression Region	1-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.