



# Recombinant Mouse Cysteine protease ATG4D (Atg4d)

<b>Product Code</b>	CSB-MP804331MO
<b>Abbreviation</b>	Atg4d
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8BGV9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNSVSPAAQ YRSGSSEDAR RADCRRPRGQ TRIPDPSNLG PSGSGVAALG SSGTDPAEPD EVDKFKAKFL TAWNNVKYGW AVKSRTSFSK ISTVHLCGRC YHFEGEGDIQ RFQRDFVSRL WLTYRRDFPP LAGGSLTSDC GWGCMLRSGQ MMLAQGLLLH FLPRDWRWVE GTGLASSEMP GPASPSRCRG PGRRGPPRWT QGALEMEQDR WHRRIVSWFA DHPRAPFGLH RLVELGRSSG KKAGDWYGPS VVAHILRKAV ESCSEVSRLV VYVSQDCTVY KADVARLLSW PDPTAEWKS VILVPVRLGG ETLNPVYVPC VKELLRSELC LGIMGGKPRH SLYFIGYQDD FLLYLDPHYC QPTVDVSQPS FPLESFHCTS PRKMAFAKMD PSCTVGFYAG NRKEFETLCS ELMRILSSSS VTERYPMFTV AEGHAQDHSL DALCTQLSQP TLRLPCTGRL LKAKRPSSD FVFL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Atg4d
<b>Protein Names</b>	Recommended name: Cysteine protease ATG4D EC= 3.4.22.- Alternative name(s): AUT-like 4 cysteine endopeptidase Autophagin-4 Autophagy-related cysteine endopeptidase 4 Autophagy-related protein 4 homolog D
<b>Expression Region</b>	1-474
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	Autophagy is the process by which endogenous proteins and damaged organelles are destroyed intracellularly. Autophagy is postulated to be essential for cell homeostasis and cell remodeling during differentiation, metamorphosis, non-apoptotic cell death, and aging. Reduced levels of autophagy have been described in some malignant tumors, and a role for autophagy in controlling the unregulated cell growth linked to cancer has been proposed. This gene encodes



a member of the autophagin protein family. The encoded protein is also designated as a member of the C-54 family of cysteine proteases.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.