



# Recombinant Mouse SPARC-related modular calcium-binding protein 1 (Smoc1)

<b>Product Code</b>	CSB-BP804372MO
<b>Abbreviation</b>	Smoc1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8BLY1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	HRTTG PRFLISDRDP PCNPHCPRTQ PKPICASDGR SYESMCEYQR AKCRDPALAV VHRGRCKDAG QSKCRLERAQ ALEQAKKPQE AVFVPECGED GSFTQVQCHT YTGycwcvtp DGKPISSGSSV QNKTPVCSPG VTDKPLSQGN SGRKVSFRFF LTLNSDDGSK PTPTMETQPV FDGDEITAPT LWIKHLVIKD SKLNNTNVRN SEKVHSCDQE RQSALEEARQ NPREGIPIE CAPGGLYKPV QCHQSTGYCW CVLVDTGRPL PGTSTRYVMP SCESDARAKS VEADDPFKDR ELPGCPEGKK MEFITSL LDA LTTDMVQAIN SAAPTGGGRF SEPDPSTLE ERVAHWYFSQ LDSNSSDDIN KREMKPKRY VKKKAKPKKC ARRFTDYCDL NKDKVISLPE LKGCLGVSKE GGSLGSFPQG KRAGTNPFIG RLV
<b>Source</b>	Baculovirus
<b>Target Names</b>	Smoc1
<b>Protein Names</b>	Recommended name: SPARC-related modular calcium-binding protein 1 Alternative name(s): SPARC-related gene protein Secreted modular calcium-binding protein 1 Short name= SMOC-1
<b>Expression Region</b>	26-463
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.