



Recombinant Mouse Voltage-dependent L-type calcium channel subunit beta-2 (Cacnb2)

Product Code	CSB-BP804482MO
Abbreviation	Cacnb2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8CC27
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MVQSDTSKSP PVAAVAQESQ MELLESAAPA GALGAQSYGK GARRKNRFKG SDGSTSSDTT SNSFVRQGS A DSYTSRPSDS DVSLEEDREA VRREAERQAQ AQLEKAKTKP VAFVRTNVR YSAAQEDDVP VPGMAISFEA KDFLHVKEKF NNDWWIGRLV KEGCEIGFIP SPVKLENMRL QHEQRAKQ GK FYSSKSGGNS SSSLGDIVPS SRKSTPPSSA IDIDATGLDA EENDIPANHR SPKPSANSVT SPHSKEKRMP FFKKTEHTPP YDVPSPMRPV VLVGPSLKG Y EVT DMMQKAL FDFLKHRFEG RISITRV TAD ISLAKRSVLN NPSKHAIER SNTRSSLAEV QSEIERIFEL ARTLQLVLD ADTINHPAQL SKTSLAPIIV YVKISSPKVL QRLIKSRGKS QAKHLNVQMV AADKLAQCPP QESFDVILDE NQLEDACEHL ADYLEAYWKA THPPSGNLPN PLLSRTLASS TLPLSPTLAS NSQGSQGDQR PDRSAPRSAS QAEEPCLEP VKKSQHRSSS ATHQNHRSGT GRGLSRQETF DSETQESRDS AYVEPKEDYS HEHVD RYVPH REHNHREETH SSSNGHRHRES RHRSRDMGRD QDHNECIKQR SRHKS KDRYC DKEGEVISKR RNEAGEWNRD VYIRQ
Source	Baculovirus
Target Names	Cacnb2
Protein Names	Recommended name: Voltage-dependent L-type calcium channel subunit beta-2 Short name= CAB2 Alternative name(s): Calcium channel voltage-dependent subunit beta 2
Expression Region	1-655
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a subunit of a voltage-dependent calcium channel protein which is a member of the voltage-gated calcium channel superfamily. The gene product was originally identified as an antigen target in Lambert-Eaton



myasthenic syndrome which is an autoimmune disorder. Mutations in this gene are associated with Brugada syndrome. Alternatively spliced variants have been identified for this gene.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.