



# Recombinant *Xenopus tropicalis* Forkhead box protein A1 (foxa1)

<b>Product Code</b>	CSB-MP804751XBF
<b>Abbreviation</b>	foxa1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8AWH1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus tropicalis</i> (Western clawed frog) ( <i>Silurana tropicalis</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLGIVKMEGH ETTDWSNYYQ DTQEAYSSVP VSNMTQGLAS MNTYMTMNPMS SSSSNMTAGS FNMSYANSGL GAGLSPSGMS GMGAGSASAM NGMGSGVSSM GTALSPSSMN AMSAQQASIN SLSYSGMNPG MSPMAYGPSN MNRTRDTKTF RRSYPHAKPP YSYISLITMA IQQAPSKMLT LSEIQWIMD LFLYYRQNNQ RWQNSIRHSL SFNDCFVKVA RSPDKPGKGS YWTLHPDSGN MFENGCYLRR QKRFKCEKQQ GGKGSQDGRK DVSGPSSPLH RVHGKSSQMD SSSSMSNPSS SPQSLEHNGS NGEMKPQVAA GPSPLSSHQN HSTHSLAHET HIHLKGDPHY SFNHFPFSINN LMSSEQQHK LDFKAYEQAL QQYSSYSGGL PGMPLGSPSM AGRGSIEPSA LEPTYQQGVY SRPVLNTS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	foxa1
<b>Protein Names</b>	Recommended name: Forkhead box protein A1 Short name= FoxA1 Alternative name(s): tFoxA1
<b>Expression Region</b>	1-428
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.