



Recombinant Mouse E3 ubiquitin-protein ligase RNF187 (Rnf187)

Product Code	CSB-MP804774MO
Abbreviation	Rnf187
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BFX1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MALPAGPADA ICALCQRAPR EPVRADCGHR FCRACVVRFW AEEDGPFPCP ECADDCWQRA VEPSRPPLSR RLLALEEAAA APARDGPASE AALQLLCRAD GDPLCSACRM AAGPEPPEWE PRWRKALRGK ENKGSVEIMR KDLNDARDLH GQAESAAAVW KGHVMDRRKK ALTDYKKLRA FFVEEEEHFL QEAEKDEGAS EDDELADPAD RFRSLLQAVS ELEKKHRNLG LSMLLQ
Source	Mammalian cell
Target Names	Rnf187
Protein Names	Recommended name: E3 ubiquitin-protein ligase RNF187 EC= 6.3.2.- Alternative name(s): RING domain AP1 coactivator 1 Short name= RACO-1 RING finger protein 187
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.