



# Recombinant Mouse RNA binding protein fox-1 homolog 2 (Rbfox2)

<b>Product Code</b>	CSB-MP804837MO
<b>Abbreviation</b>	Rbfox2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8BP71
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAEGGQAQQQ PPQLGPGAAA RGMKRESEVE LPVPGAGADG PEPGLSKRPR TEEAADGGMQ NEPLTPGYHG FPARDGQGNQ EPTTTPDAMV QPFTTIPFPP PPQNGIPTEY GVPHTQDYAG QTSEHNLTLY GSTQPHGEQS SNSPSNQNGS LTQTEGGAQT DGQQSQTQSS ENSEKSTPK RLHVSNIPIFR FRDPDLRQMF GQFGKILDVE IIFNERGSKG FGVTFENSA DADRAREKLH GTVVEGRKIE VNNATARVMT NKKMVTPIYAN GWKLSVVGA VYGPELYAAS SFQADVSLGN EAAVPLSGRG GINTYIPLIS LPLVPGFPYP TAATTAAAFR GAHLRGRGRT VYGAVRAVPP TAIPAYPGVV YQDGFYGADL YGGYAAARYA QPATATAATA AAAAAAYSD GYGRVYTADP YHALAPAASY GVGAVASLYR GGYSRFAPY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Rbfox2
<b>Protein Names</b>	Recommended name: RNA binding protein fox-1 homolog 2 Alternative name(s): Fox-1 homolog B Fox-1 homolog Fxh Hexaribonucleotide-binding protein 2 RNA-binding motif protein 9 RNA-binding protein 9
<b>Expression Region</b>	1-449
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.