



Recombinant Mouse RNA binding protein fox-1 homolog 2 (Rbfox2)

Product Code	CSB-YP804837MO
Abbreviation	Rbfox2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BP71
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MAEGGQAQQQ PPQLGPGAAA RGMKRESEVE LPVPGAGADG PEPGLSKRPR TEEAADGGMQ NEPLTPGYHG FPARDGQGNQ EPTTTPDAMV QPFTTIPFPP PPQNGIPTEY GVPHTQDYAG QTSEHNLTLY GSTQPHGEQS SNSPSNQNGS LTQTEGGAQT DGQQSQTQSS ENSEKSTPK RLHVSNIPIFR FRDPDLRQMF GQFGKILDVE IIFNERGSKG FGVTFENSA DADRAREKLH GTVVEGRKIE VNNATARVMT NKKMVTPIYAN GWKLSVVGA VYGPELYAAS SFQADVSLGN EAAVPLSGRG GINTYIPLIS LPLVPGFPYP TAATTAAAFR GAHLRGRGRT VYGAVRAVPP TAIPAYPGVV YQDGFYGADL YGGYAAARYA QPATATAATA AAAAAAYSD GYGRVYTADP YHALAPAASY GVGAVASLYR GGYSRFAPY
Source	Yeast
Target Names	Rbfox2
Protein Names	Recommended name: RNA binding protein fox-1 homolog 2 Alternative name(s): Fox-1 homolog B Fox-1 homolog Fxh Hexaribonucleotide-binding protein 2 RNA-binding motif protein 9 RNA-binding protein 9
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.