



# Recombinant Mouse Rab11 family-interacting protein 4 (Rab11fip4)

<b>Product Code</b>	CSB-BP804842MO
<b>Abbreviation</b>	Rab11fip4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8BQP8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAGGAGWAGA PAALLRSVRR LREVFVCGR DPDGFLRVER VAALGLRFGQ GEEVEKLVKC LDPNDLGRIN FKDFCRGVFA MKGCEELLKD VLSVESAGTL PCSPDIPDCV EQGSDFSGST DGEQLPREPD FFQEDEEEAM TLALPEGPQE LDMDSPMESS QGPEGSVKGC GEEKEPELGG LFLPEDKCLV LTPSVTTSDL STHSTASLIS NEEQFEDYGE GDDVDCAPSS PCPDDETR TN VYSDLGSSVS SSAGQTPRKM RHAYNSELLD VYCSQCCKKI NLLNDLEARL KNLKANS PNR KISSTAFGRQ LMHNSNFSSS NGSTEDLFRD SIDSCDNDIT EKVSFLEKKV TELENDLSTS GGLKSKLKQE NMQLVHRVHE LEEMVKDQET TAEQALEEEA RRHREVYCKL EREKSTELEL LNTRVQQLEE ENTDLRRTVA RLKSQTEKLD EERQRMSDRL EDTSLRLKDE MDLYKRMMDK LRQNRLEFQK EREA TQELIE DLRRELEHLQ MYKLCDCERPG RGRSSSGLGE FNARAREVEL EHEVKRLKQE NHKLRDQNDL LNGQILSLSL YEAKNLFATQ TKAQSLAAEI DTASRDELME ALKEQEEINF RLRQYMDKII LAILDHNPSI LEIKH
<b>Source</b>	Baculovirus
<b>Target Names</b>	Rab11fip4
<b>Protein Names</b>	Recommended name: Rab11 family-interacting protein 4 Short name= FIP4- Rab11 Short name= Rab11-FIP4 Short name= mRab11-FIP4A
<b>Expression Region</b>	1-635
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.