



Recombinant Mouse Rab11 family-interacting protein 4 (Rab11fip4)

Product Code	CSB-EP804842MO
Abbreviation	Rab11fip4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8BQP8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MAGGAGWAGA PAALLRSVRR LREVFVCGR DPDGFLRVER VAALGLRFGQ GEEVEKLVKC LDPNDLGRIN FKDFCRGVFA MKGCEELLKD VLSVESAGTL PCSPDIPDCV EQGSDFSGST DGEQLPREPD FFQEDEEEAM TLALPEGPQE LDMDSPMESS QGPEGSVKGC GEEKEPELGG LFLPEDKCLV LTPSVTTSDL STHSTASLIS NEEQFEDYGE GDDVDCAPSS PCPDDETR TN VYSDLGSSVS SSAGQTPRKM RHAYNSELLD VYCSQCCKKI NLLNDLEARL KNLKANS PNR KISSTAFGRQ LMHNSNFSSS NGSTEDLFRD SIDSCDNDIT EKVSFLEKKV TELENDLTS GGLKSKLKQE NMQLVHRVHE LEEMVKDQET TAEQALEEEA RRHREVYCKL EREKSTELEL LNTRVQQLEE ENTDLRRTVA RLKSQTEKLD EERQRMSDRL EDTSLRLKDE MDLYKRMMDK LRQNRLEFQK EREAQELIE DLRRELEHLQ MYKLCDCERPG RGRSSSGLGE FNARAREVEL EHEVKRLKQE NHKLRDQNDL LNGQILSLSL YEAKNLFATQ TKAQSLAAEI DTASRDELME ALKEQEEINF RLRQYMDKII LAILDHNPSI LEIKH
Source	E.coli
Target Names	Rab11fip4
Protein Names	Recommended name: Rab11 family-interacting protein 4 Short name= FIP4- Rab11 Short name= Rab11-FIP4 Short name= mRab11-FIP4A
Expression Region	1-635
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.