



Recombinant Mouse WD repeat-containing and planar cell polarity effector protein fritz homolog (Wdpcp)

Product Code	CSB-EP804911MO-B
Abbreviation	Wdpcp
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8C456
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSFCLTELHL WSLKSTLHIA DRDIGVYQYY DKKDPSVSAT EHGNLEEKQR LAESRDYPWT LKNRRPEKLR DSLKELEELM QSSPCVLSKW KSKYICQLLF GSGVLVLSL SGPQLEKVI DRSLVGKLIS DTISDALLTD SFIILSFLAQ NKLCFIQFTK KMDSL DGNKR LEKLSALDLK ISYYDIPGPA NRTIDRHLAV NSTQDLVVCW WPLVSDDVWP WTPVSSEKDR ANMLLLGFTQ GGLEVLSFVR TEWSPLDVHF GTKQPYQVFT VECSVSVDKE PMADSCIYES VRNKLHCVSV TRIPLRSKAI SCCRNSTEDK LIVGCEDSSV ILYEAHRGVT LLAQAE LRPS LISCHPSGAI LLVGSNQGEL QIFDIALSPI NIQLLAEDYS PKETLQFKKF FDVSSSLVQM QWMAPPVVFQ KPKRGEICDL LFLRFNKGPL GVLLFKLGIL TRGQLGLVDL ILQYIHYSEV YE AISILRSM DWDTLGQQCL IGMG TIVNHL LRQRLTPERE AQLEASLGTF YAPTRPLD T TLEYREPVS KYARRLFHHL LRYKRFEKAF LLAVDIGARD LFMDIHYLAL DMGELALAEV ARRRAHDIDV ESVSSGV ELL GPLDRRDMLN EGFASSALMP EGENKFPGLL PSIGSTHMQT LQQKIPNGPS SRWAIERTE EEEEEEEEEEE EELCTDSSGA TTWNAEGELK EDQRKQDIGD VGSLRMVHFG LV
Source	E.coli
Target Names	Wdpcp
Protein Names	Recommended name: WD repeat-containing and planar cell polarity effector protein fritz homolog Short name= mFrtz Alternative name(s): Homoloc-13 WD repeat-containing and planar cell polarity effector protein
Expression Region	1-722
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.