



Recombinant Mouse WD repeat-containing protein 24 (Wdr24)

| | |
|--------------------------|--|
| Product Code | CSB-BP804957MO |
| Abbreviation | Wdr24 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8CFJ9 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MEKMSRVSTA LSGSALTGRT MHCHLDAPAN AISVCRDAAQ VVAGRSIFK IYAIEEEQFV EKLNLRVGRK PSLNLSCADV VWHQMDENLL ATAATNGVVV TWNLGRPSRN KQDQLFTEHK RTVNVKCFHP TEAHVLLSGS QDGFMKCFDL RRKDSVSTFS GQSESVRDVQ FSIRDYFTFA STFENGNVQL WDIRRPDRCE RMFTAHNQPV FCCDWHPEDR GWLATGGRDK MVKVWDMTTH RAKEIHCVQT IASVARVKWR PECRHHLATC SMMVDHNIYV WDVRRPFVPA AMFEEHRDVT TGIWRPHPHD PSFLLSGSKD STLCQHFRD ASQPVERANP EGLCYGLFGD LAFVAVKESLV AESGRKPYA GDRRHPIFFK RKLDPAEFVS GLASSALSVF ETESSGGSMS WFVDTAERYV LAGRPLAELC DHNAKVAREL GRNQVAQTWT MLRIYCSPG LVSSANLNHS VGKGSSCGLP LMNSFNLKDM GPGLGSETRL DRSKGDTRSD AALLDSSATL VTNEDNEETE GSDVPADYLL GDVEGEDDEL YPLDTEHVHS EEPEYVLPQE AFPLRHEIVD TPSGPEHLQD KADSPHVSGN EADTASLAPV DSSSSLLSVS HALYDSRLPP DFFSVLVRDM LRFYAEQGDV QMAVSVLIVL GERVRKDIDE QTQEHWYTSY IDLLQRFCLW NVSNEVVKLS TSRAVSLNQ ASTTLHVNCS HCKRPMSSRG WVCDRCHRCA SMCAVCHHVV KGLFVWCQGC SHGGHLQHIM KWLEGSSHCP AGCGHLCEYS |
| Source | Baculovirus |
| Target Names | Wdr24 |
| Protein Names | Recommended name: WD repeat-containing protein 24 |
| Expression Region | 1-790 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a |



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.